

UNLOCKING OPTIMAL CONDITIONS FOR HEK293T CELL EXPANSION USING SCINUS BIOREACTOR TECHNOLOGY

Javier Olmos Becerra^{1,2} (javier.olmos@scinus.com), Iris Pijnenburg¹ (iris.pijnenburg@scinus.com), Orsolya Frittmann² (o.frittmann@umcg.nl), Edwin Bremer² (e.bremer@umcg.nl) and Ruud Das¹ (ruud.das@scinus.com)

1. Scinus Cell Expansion Netherlands B.V., Bilthoven, Netherlands

2. Department of Hematology, University of Groningen, University Medical Center Groningen (UMCG), Groningen, Netherland



INTRODUCTION

HEK293T cells are a versatile cell line that can, among others, be used as a producer cell for lentiviral vectors (LV) for cell and gene therapy. In particular, HEK293T can be used for LV production for Chimeric Antigen-Receptor T cell (CAR-T) therapy. Innovations in CAR-T therapy might include those that facilitate point-of-care production, instead of logistically challenging centralized production. Microcarriers are an appealing culture substrate due to their efficient surface-to-volume ratio, but efficient HEK293T culture in agitation systems is challenging because of their sensitivity to shear stresses. In this study we explore the ability to employ microcarrier-based HEK293T cultivation for use in closed bioreactor technology (SCINUS bioreactor) for point-of-care LV production.

METHODS AND MATERIALS

Design-of-experiments for HEK293 culture on Cytodex-1

A bag-in-bag system was developed, to test multiple conditions with the agitation approach of the SCINUS system. A two-level factorial design was created for three parameters: rocking speed, seeding density and static pause. PDT (day 1-3) was used as a response.

	Low	Center	High
Rocking speed (deg/sec)	60	190	320
Seeding density (cells/mL)	5×10^4	1.25×10^5	2×10^5
Static interval (H)	1	2.5	4

Expansion and bead to bead transfer

The ability of HEK293 cells to proliferate under agitated conditions and subsequently display bead-to-bead transfer was investigated using two different microcarrier cultures in spinners. Cells were seeded at various seeding densities (50 – 200k cells/mL). Volume and microcarrier content were increased 9-fold on day 3. Growth and bead-to-bead transfer were then assessed on day 7.

SCINUS culture and transfectability test

HEK293T cells were cultured in the SCINUS bioreactor using LDC-DMC. After 8 days of culture, cells were harvested and then seeded for transfection. Transfectability of the cells was evaluated by using a GFP plasmid and then examining the expression through flow cytometry.

RESULTS

DoE experiment

No significant interaction was found between the studied parameters. Response surface models (Figure 2) indicated that a lower seeding density ($p=0.0002$) and a lower initial static pause ($p=0.0162$) resulted in faster cell growth. Distribution was comparable with the spinner cultures (Figure 3)





Figure 2 – Response surface model for the three different rocking speeds.

SPINNER



SCINUS




Figure 3 – Visual inspections at day 3 of HEK293T cells on Cytodex 1 (3 g/L) cultured at different seeding densities using spinners and the SCINUS bioreactor.

Bead-to-bead transfer

Both microcarriers supported the growth of HEK293T cells, but Cytodex1 showed clumped growth on old microcarriers, while in LDC-DMC cells were well distributed, showing bead-to-bead transfer (Figure 4).

CYTODEX 1



LDC-DMC





Figure 4 – Visual inspections at day 7 (harvest day) of HEK293T-17 on Cytodex1 and LDC-DMC.

Expansion in SCINUS and transfection

HEK293T-17 cells demonstrated the ability to proliferate in the SCINUS bioreactor, obtaining >800 million cells in 8 days. Aggregation was observed (Figure 5A) but harvested and transfected cells showed high level of GFP expression (Figure 5B), comparable to the monolayer control ($p=0.089$).

(A)



(B) Transfectability after 48h




Figure 5 – (A) Visual inspections at day 8 (harvest day). (B) GFP expression levels after 48h.

CONCLUDING REMARKS

This study highlights the integration of microcarrier-based HEK293T culture into the SCINUS bioreactor, laying the groundwork for point-of-care lentiviral vector production for advanced therapies, such as CAR-T cell therapy.




Figure 1 – General scheme of the experiments.